

Product datasheet

Recombinant Human DCUN1D4 protein ab171715

1 Image

Description

Product name	Recombinant Human DCUN1D4 protein
Purity	> 85 % SDS-PAGE. ab171715 was purified by anion-exchange and gel-filtration chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q92564</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MGSSHHHHHH SSGLVPRGSH MGSMHSDAAA VNFQLNSHLS TLANIHKIYH TLNKLNLTED IGQDDHQTGS LRSCSSSDCF NKVMPPRKKR RPASGDDLSA KKSRRHDSMYR KYDSTRIKTE EAAFSSKRCL EWFYEYAGTD DVVGPEGMEK FCEDIGVEPE NVVMLVLAWK LDAQNMGYFT LQEWLKGMS LQCDTTEKLR NTLDYLRSL NDSTNFKLIY RYAFDFAREK DQRS LDINTA KCMLG LLLGK WPLFPVFHQ FLEQSKYKVI NKDQWCNVLE FSRTINLDLS NYDEDGAWPV LLDEFVEWYK DKQMS</p>
Predicted molecular weight	37 kDa including tags
Amino acids	1 to 292
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab171715** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8


Constituents: 0.03% DTT, 0.32% Tris HCl, 0.06% EDTA, 50% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

General Info

Sequence similarities

Contains 1 DCUN1 domain.

Images



15% SDS-PAGE analysis of ab171715 (3µg).

SDS-PAGE - Recombinant Human DCUN1D4
protein (ab171715)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors