

Recombinant Human DPF2/REQ protein ab159327

1 Image

Description

Product name	Recombinant Human DPF2/REQ protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MAAVVENVVKLLGEQYYKDAMEQCHNYNARLCAERSVRL PFLDSQTGVAQ SNCYWMEKRHRGPGLASGQLYSYPARRWRKKRAHPP EDPRLSFPSIKP DTDQTLKKEGLISQDGSSLEALLRTDPLEKRGAPDPRVD DDSLGEFPVTN SRARKRILEPDDFLDDLDDDEYEDTPKRRGKGKSKGKG VGSARKKLDAS ILED RD KPYACDICGKRYKNRPGLSYHYAHSHLAE E EGED KEDSQPPTPV SQRSEEQSKKGPDGLALPNNYCDFCLGDSKINKKTGQP EELVSCSDCGR SGHP SCLQFTPVMMAAVKTYRWQCIECKCCNICGTSEND DQLLFCDDCDR GYHMYCLTPSMSEPPEG SWSCHLCLDLLKEKASYQNQN SS
Amino acids	1 to 391
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab159327** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid

Additional notes

This product was previously labelled as DPF2.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

May be a transcription factor required for the apoptosis response following survival factor withdrawal from myeloid cells. Might also have a role in the development and maturation of lymphoid cells.

Tissue specificity

Ubiquitous.

Sequence similarities

Belongs to the requiem/DPF family.


Contains 1 C2H2-type zinc finger.

Contains 2 PHD-type zinc fingers.

Cellular localization

Nucleus. Cytoplasm. 30% nuclear. 70% cytoplasmic.

Images



SDS-PAGE - Recombinant Human DPF2/REQ protein (ab159327)

ab159327 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors