

Recombinant Human FAM84B protein ab174398

1 Image

Description

Product name	Recombinant Human FAM84B protein
Purity	> 85 % SDS-PAGE. ab174398 is purified using conventional chromatography.
Expression system	Escherichia coli
Accession	<u>Q96KN1</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSMGNQVEKLTHLSYKE VPTADPTGVDRD DGPRIGVSYIFSNDDEDVEPQPPQGPDGGGLPDGGDG PPPPQPQPYDPR LHEVECSVFYRDECIYQKSFAPGSAALSTYTPENLLNKCK PGDLVEFVSQ AQYPHWAVYVGNFQVVHLHRLEVINSFLTDASQGRRGRV VNDLYRYKPLS SSAVVRNALAHVGAKERELSWRNSESFAAWCRYGKREF KIGGELRIGKQP YRLQIQLSAQRSHTLEFQSLEDLIMEKRRNDQIGRAAVLQE LATHLHPAE PEEGDSNVARTTPPPGRPPAPSSEEDGEAVAH
Predicted molecular weight	37 kDa including tags
Amino acids	1 to 310
Tags	His tag N-Terminus
Additional sequence information	NCBI Accession No. NP_777571.

Specifications

Our **Abpromise guarantee** covers the use of **ab174398** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry
--------------	-------------------

	SDS-PAGE
<b>Mass spectrometry</b>	MALDI-TOF
<b>Form</b>	Liquid


## Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.02% DTT, 0.002% PMSF, 0.88% Sodium chloride
------------------------------	---

## General Info

<b>Tissue specificity</b>	Expressed in esophageal squamous cell carcinomas.
<b>Sequence similarities</b>	Belongs to the FAM84 family.

## Images



15% SDS-PAGE analysis of ab174398 (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors