

Recombinant Human FIP protein ab152790

1 Image

Description

Product name	Recombinant Human FIP protein
Expression system	Wheat germ
Accession	<u>Q15853</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MDMLDPGLDPAASATAAAAAASHDKGPEAEEGVELQEGG DGPGAEEQTAVA ITSVQQAAGFDHNIQYQFRTE TNGGQVTYRVVQVTDGQLD GQGDTAGAVS VVSTAAGAGGQQAVTQVGVDGAAQRPGPAAASVPPGPA APFPLAVIQNPF SNGGSPAAEAVSGEARFAYFPASSVGDTTAVSVQTTDQS LQAGGQFYVMM TPQDVLQTGTQRTIAPRTHPYSPKIDGTRTPRDERRRAQH NEVERRRRDK INNWIQLSKIIPDCNADNSKTGASKGGILSKACDYIRELRQ TNQRMQET FKEAERLQMDNELLRQQIEELKNENALLRAQLQQHNLEM VGEGTRQ
Predicted molecular weight	63 kDa including tags
Amino acids	1 to 346

Specifications

Our <b>Abpromise guarantee</b> covers the use of <b>ab152790</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	ELISA SDS-PAGE Western blot
Form	Liquid

## Additional notes

Protein previously labeled as USF2.

## Preparation and Storage

### Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info

### Function

Transcription factor that binds to a symmetrical DNA sequence (E-boxes) (5'-CACGTG-3') that is found in a variety of viral and cellular promoters.

### Tissue specificity

Ubiquitous.


### Sequence similarities

Contains 1 basic helix-loop-helix (bHLH) domain.

### Cellular localization

Nucleus.

## Images



12.5% SDS-PAGE showing ab152790 stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors