

Product datasheet

Recombinant Human FLJ25006 protein ab203554

1 Image

Description	
Product name	Recombinant Human FLJ25006 protein
Purity	> 70 % Densitometry. Affinity purified.
Expression system	Baculovirus infected Sf9 cells
Accession	<u>Q96LW2</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGAVSCRQGQHTQQGEHTRVAVPHKQGGNIRGPWARG WKS LWTGLGTIRS DLEELWELRGHHYH QESLKPAPVLVEKPLPEWPVPQFI NLFLPEFP IRP IRGQQQLKILGLVAKGSFGTVLKVLDCTQKAVFAVKVVPK VKVLQRD TVR QCKEEVSIQRQINHPFVHSLGDSWQGKRHL FIMCSYCSTDLYSLWSAV GCFPEASIRLFAAELVLVLCYLHDLGIMHRDVKMENILLDE RGHLKL TDF GLSRHVPQGAQAYTICGTLQYMGERG
Predicted molecular weight	57 kDa including tags
Amino acids	1 to 274
Tags	GST tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab203554 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped on Dry Ice. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.50


Constituents: 0.79% Tris HCl, 0.87% Sodium chloride, 0.31% Glutathione, 0.003% EDTA, 0.0038% DTT, 0.00174% PMSF, 25% Glycerol (glycerin, glycerine)

General Info

Relevance

FLJ25006 likely exerts a protein serine/threonine kinase activity.

Images



SDS-PAGE - Recombinant Human FLJ25006 protein
(ab203554)

Analysing the purity of ab203554 by densitometry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors