

Product datasheet

Recombinant Human GCDFP 15 protein (His tag)
ab235867

1 Image

Description

Product name	Recombinant Human GCDFP 15 protein (His tag)		
Purity	> 90 % SDS-PAGE.		
Expression system	Escherichia coli		
Accession	<u>P12273</u>		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	QDNTRKIIIIKNFDIPKSVRPNDEVTAVLAVQTELKECMVVKT YLISSIPL QGAFNYKYTACLCDDNPKTFYWDFYTNRTVQIAAVVDVIR ELGICPDDAA VIPIKNNRFYTIEILKVE		
Predicted molecular weight	18 kDa including tags		
Amino acids	29 to 146		
Tags	His tag N-Terminus		
Additional sequence information	Full-length mature chain lacking the signal peptide. N-terminal 6xHis-tagged.		

Specifications

Our <u>Abpromise guarantee</u> covers the use of ab235867 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.2
-----------------------	--

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

General Info

Tissue specificity	Expressed in pathological conditions of the mammary gland and in several exocrine tissues, such as the lacrimal, salivary, and sweat glands.
Sequence similarities	Belongs to the PIP family.
Cellular localization	Secreted.

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel of ab235867.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors