

Product datasheet

Recombinant Human HHLA2 protein (Tagged) ab271384

[1 Image](#)

Description

Product name	Recombinant Human HHLA2 protein (Tagged)
Purity	>= 90 % SDS-PAGE.
Expression system	HEK 293 cells
Accession	<u>Q9UM44</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	IFPLAFFIVPMNEQVIGRLDEDIILPSSFERGSEVVIHWKY QDSYKVH SYYKG SDHLESQDPRYANRTSLFYNEIQNGNASLFFRRVSLLDGI YT CYVGTAIQVITN KVVLLKVGVFLTPVMKYEKRNTNSFLICSVLSVYPRP IITWKMDNTPISENNMEET GSLDSFSINSPLNITGSNSSYECTIENSL LKQWTWTGRWTMKDGLHKMQSEHVLSL CQPVNDYFSPNQDFKVTWSRMK SGTFSVLAYLSSSQNTIINESRFSWNKELINQ SDFSMNLMDLNLSDS GEYLCNISSDEYLLTIHTVHVEPSQETASHNIEGRMDPK SCDKTHTC PPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDV SHEDPEV K FNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDW LNGKEYKCKVSNKAL PAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLV KGF YPSDIAVEWES NGQPENNYKTTPPVLDSDGSFFLYSKLTVDKSRWQQG NVFSCSVMHEALHNHYTQ KSLSLSPGKGGGLNDIFEAQKIEWHE
Predicted molecular weight	66 kDa

Amino acids 23 to 344
Tags Avi tag C-Terminus , Fc tag C-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab271384** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE
Form Liquid
Additional notes Enzymatically biotin-labeled using Avi-tag™ technology


Preparation and Storage

Stability and Storage Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.
pH: 7.40
Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

General Info

Tissue specificity Expressed in colon, small intestine, lung and kidney.
Sequence similarities Contains 1 Ig-like C1-type (immunoglobulin-like) domain.
Contains 2 Ig-like V-type (immunoglobulin-like) domains.
Cellular localization Membrane.

Images



SDS-PAGE analysis of ab271384.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors