

Product datasheet

Recombinant Human Homeobox protein SIX6
ab134617

1 Image

Description	
Product name	Recombinant Human Homeobox protein SIX6
Purity	> 90 % SDS-PAGE. ab134617 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>O95475</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSMFQLPILNFSPQQVA GVCETLEESGDV ERLGRFLWSLPVAPAACEALNKNESVLRARAIVAFHGGN YRELYHILENH KFTKESHAKLQALWLEAHYQEA EKLRGRPLGPVDKYRVR KKFPLPRTWD GEQKTHCFKERTRNLLREWYLQDPYPNPSKKRELAQATG LTPTQVGNWFK NRRQRDRAAAAKNRLQQQVLSQGSGRALRAEGDGTPEV LGVATSPAASLS SKAATSAISITSSDSECDI
Predicted molecular weight	30 kDa including tags
Amino acids	1 to 246
Tags	His tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab134617 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF

Form Liquid

Preparation and Storage

Stability and Storage


Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.08% DTT, 0.32% Tris HCl, 0.06% EDTA, 40% Glycerol (glycerin, glycerine), 2.34% Sodium chloride

General Info

Images



15% SDS-PAGE analysis of 3 µg ab134617.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors