

Recombinant Human IL-17B protein ab50161

Description

| | |
|--------------------------|---|
| Product name | Recombinant Human IL-17B protein |
| Purity | > 98 % SDS-PAGE. Greater than 98% by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g). |
| Expression system | Escherichia coli |
| Protein length | Full length protein |
| Animal free | No |
| Nature | Recombinant |
| Species | Human |
| Sequence | MQPRSPKSKR KGQGRPGPLA PGPHQVPLDL VSRMKPYARM EEYERNIEEM VAQLRNSSEL AQRKCEVNLQ LWMSNKRSLS PWGYSINHDP SRIPVDLPEA RCLCLGCVNP FTMQEDRSMV SVPVFSQVPV RRRLCPPPPR TGPCRQRAVM ETIAVGCTC I F |

Specifications

Our **Abpromise guarantee** covers the use of **ab50161** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|---------------------|--------------------------------|
| Applications | SDS-PAGE Functional Studies |
| Form | Lyophilized |

Preparation and Storage

| | |
|------------------------------|--|
| Stability and Storage | Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term. |
| Reconstitution | For lot specific reconstitution information please contact our Scientific Support Team. |

General Info

| | |
|------------------------------|---|
| Function | Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1. |
| Tissue specificity | Expressed in adult pancreas, small intestine, stomach, spinal cord and testis. Less pronounced expression in prostate, colon mucosal lining, and ovary. |
| Sequence similarities | Belongs to the IL-17 family. |
| Cellular localization | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors