

## Product datasheet

# Recombinant Human ISLR protein ab127069

### Description

<b>Product name</b>	Recombinant Human ISLR protein
<b>Purity</b>	> 90 % SDS-PAGE. Purified via His tag
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<b><u>O14498</u></b>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Predicted molecular weight</b>	18 kDa
<b>Amino acids</b>	111 to 275

### Specifications

Our **Abpromise guarantee** covers the use of **ab127069** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE
<b>Form</b>	Lyophilized

### Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Store at -20°C.  Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
<b>Reconstitution</b>	Reconstitute with water to desired concentration.

### General Info

<b>Tissue specificity</b>	Expressed in various tissues including retina, heart, skeletal muscle, prostate, ovary, small intestine, thyroid, adrenal cortex, testis, stomach and spinal cord.
<b>Sequence similarities</b>	Contains 1 Ig-like (immunoglobulin-like) domain. Contains 5 LRR (leucine-rich) repeats.

	Contains 1 LRRCT domain.
	Contains 1 LRRNT domain.
<b>Developmental stage</b>	Expressed in fetal lung and kidney.
<b>Cellular localization</b>	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors