

Recombinant Human KIAA0020 protein ab160721

1 Image

Description

Product name	Recombinant Human KIAA0020 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MEVKGKKQFTGKNTKTAQEKNRFHKNSDSGSSKTFPTRK VAKEGGPKVTS RNFEKSITKLGGKGVKQFKNKQQGDKSPKNKFQPANKFN KKRKFPDGRS DESAAKPKWDDFKKKKKELKQSRQLSDKTNYDIVRAK QMWEILRRKDC DKEKRVKLMSDLQKLIQGKIKTIAFAHDSTRVIQCYIQYGNE EQRKQAFE ELRDDLVELSKAKYSRNIVKKFLMYGSKPQAEIIRSFKGHV RKMLRHAE ASAVEYAYNDKAILEQRNMLTEELYGNTFQLYKSADHPTL DKVLELQPE KLELIMDEMKQILTPMAQKEAVIKHSLVHKVFLDFFTYAPP KLRSEMIEA IREAVVYLAHTHDGARVAMHCLWHGTPKDRKVIVKTMKTY VEKVANGQYS HLVLLAAFDCIDDTKLVKQIISEISSLPSIVNDKYGRKVLLYL LSPRD PAHTVREIIEVLQKGDGNAHKKDTEVRRRELLESISPALL SYLQEHAQE VVLDKSACVLVSDILGSATGDVQPTMNAIASLAATGLHPG GKDGEIHIAE HPAGHLVLKWLEQDKKMKENGREGCFAKTLVEHVGMMKN LKSWASVNRGA IILSSLLQSCDLEVANKVKAALKSLIPTLEKTKSTSKGIEILLE KLST
Amino acids	1 to 648
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab160721** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot
Form	Liquid
Additional notes	


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Tissue specificity	Widely expressed.
Sequence similarities	Contains 1 PUM-HD domain.
	Contains 6 pumilio repeats.
Domain	The HA-8 region can be cleaved and exposed at the cell surface where it plays a role as a minor histocompatibility HLA-A*0201-restricted antigen.
Cellular localization	Nucleus > nucleolus.

Images



ab160721 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors