

Recombinant Human LRG1/LRG protein (denatured)
ab171596

1 Image

Description	
Product name	Recombinant Human LRG1/LRG protein (denatured)
Purity	> 90 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>P02750</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSVTLSPKDCQVFRSDH GSSISCQPPAEI PGYLPADTVHLAVEFFNLTHLPANLLQGASKLQELHLSSN GLESLSPEFL RPVPQLRVLDLTRNALTGLPPGLFQASATLDTLVLENQL EVLEVSWLHG LKALGHLDLSGNRLRKLPGLLANFTLLRDLGENQLETL PPDLLRGPL QLERLHLEGNKLQVLGKDLLLPQPDRLRYLFLNGNKLARVA AGAFQGLRQL DMLDLSNNSLASVPEGLWASLGQPNWDMRDGFDISGNP WICDQNLSDLYR WLQAQKDKMFSQNDTRCAGPEAVKGQTLAVAKSQ
Predicted molecular weight	36 kDa including tags
Amino acids	36 to 347
Tags	His tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab171596 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE

Form	Liquid
Additional notes	This product was previously labelled as LRG1


Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 2.4% Urea, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)
------------------------------	--

General Info

Tissue specificity	Plasma.
Sequence similarities	Contains 8 LRR (leucine-rich) repeats. Contains 1 LRRCT domain.
Cellular localization	Secreted.

Images



15% SDS-PAGE analysis of 3µg ab171596.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors