

Product datasheet

Recombinant Human Lymphotactin/ATAC protein ab85406

[1 Image](#)

Description

Product name	Recombinant Human Lymphotactin/ATAC protein
Purity	> 95 % SDS-PAGE. ab85406 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MVGSEVSDKR TCVSLTTQRL PVSRIKTYTI TEGSLRAVIF ITKRGLKVCA DPQATWVRDV VRSMDRKSNT RNNMIQTKPT GTQQSTNTAV TLTG

Specifications

Our **Abpromise guarantee** covers the use of **ab85406** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid
Additional notes	This product was previously labelled as Lymphotactin

Preparation and Storage


Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.0308% DTT, 0.316% Tris HCl, 30% Glycerol (glycerin, glycerine), 1.16% Sodium
------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

chloride

General Info

Function	Chemotactic activity for lymphocytes but not for monocytes or neutrophils.
Tissue specificity	Highest level in spleen, lower in peripheral leukocytes and very low levels in lung, colon and small intestine.
Sequence similarities	Belongs to the intercrine gamma family.
Cellular localization	Secreted.

Images



15% SDS-PAGE showing ab85406 (3µg).

Note: Real molecular weight on SDS-PAGE is shifted up compared to predicted molecular weight.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors