

Recombinant Human MA2 protein (Tagged) ab235855

1 Image

Description	
Product name	Recombinant Human MA2 protein (Tagged)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>Q9UL42</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	ALALLEDWCRIMSVDEQKSLMVTGIPADFEEAEIQEVLQE TLKSLGRYRL LGKIFRKQENANAVLLELLEDTDVSAIPSEVQGKGGVWKV IFKTPNQDTE FLERLNLFLEKEGQTVSGMFRALGQEGVSPATVPCISPEL LAHLLGQAMA HAPQPLLPMRYRKLRVFGSAVPAPEEESFEVWLEQATE VKEWPVTEAE KKRWLAESLRGPALDLMHIVQADNPSSISVEECLEAFKQVF GSLESRRTAQ VRYLKTYQEEGEKVSAYVLRLETLLRRAVEKRAIPRRIADQ VRLEQVMAG ATLNQMLWCRLRELKDQGPPPSFLELMKVIREEEEEEEAS FENESIEEPEE RDGYGRWNHE
Predicted molecular weight	57 kDa including tags
Amino acids	2 to 361
Tags	His tag N-Terminus
Additional sequence information	N-terminal 6xHis-SUMO-tagged.

Specifications

Our **Abpromise guarantee** covers the use of **ab235855** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.2

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)


General Info

Tissue specificity Brain specific. In some patients suffering from cancers, it is also specifically expressed by the testicular tumor cells.

Sequence similarities Belongs to the PNMA family.

Cellular localization Nucleus > nucleolus.

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis with 5% enrichment gel and 15% separation gel of ab235855.

SDS-PAGE - Recombinant Human MA2 protein
(Tagged) (ab235855)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors