

Recombinant Human MESDC1 protein - BSA and Azide free ab180334

1 Image

Description	
Product name	Recombinant Human MESDC1 protein - BSA and Azide free
Purity	> 95 % SDS-PAGE. ab180334 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q9H1K6</u>
Protein length	Full length protein
Animal free	No
Carrier free	Yes
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMASGSAG KPTGEAASPA PASAIGGASS QPRKRLVSVC DHCKGKMQLV ADLLLLSSEA RPVLFEGPAS SGAGAESFEQ CRDTIARTK GLSILTHDVQ SQLNMGRFGE AGDSLVELGD LVVSLTECSA HAAYLAAVAT PGAQPAQPGL VDRYRVTRCR HEVEQGCAVL RATPLADMTP QLLLEVSQGL SRNLKFLTDA CALASDKSRD RFSREQFKLG VKCMSTSASA LLACVREVKV APSELARSRC ALFSGPLVQA VSALVGFAVE PQFLGAAAV SAEGKAVQTA ILGGAMSVVS ACVLLTQCLR DLAQHPDGGA KMSDHRERLR NSACAVSEGC TLLSQALRER SSPRTLPPVN SNSVN
Predicted molecular weight	40 kDa including tags
Amino acids	1 to 362
Tags	His tag N-Terminus
Additional sequence information	(NP_072088).
Description	Recombinant Human MESDC1 protein (BSA and azide free)

Specifications

Our **Abpromise guarantee** covers the use of **ab180334** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Mass Spectrometry

SDS-PAGE

Mass spectrometry MALDI-TOF

Form Liquid

Preparation and Storage

Stability and Storage


Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 0.58% Sodium chloride, 10% Glycerol (glycerin, glycerine), 0.02% DTT

General Info

Images



15% SDS-PAGE analysis of ab180334 (3µg).

SDS-PAGE - Recombinant Human MESDC1 protein

- BSA and Azide free (ab180334)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors