

Recombinant Human MKX protein ab166113

★★★★★ 1 Abreviews 1 Image

Description	
Product name	Recombinant Human MKX protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MNTVFNKLSSAVLFEDGGASERERGGRPYSGVLDSPHA RPEVGIPDGPP LKDNLGLRHRRTGARQNGGKVRHKRQALQDMARPLKQW LYKHRDNPYPTK TEKILLALGSQMTLVQVSNWFANARRRLKNTVRQPDLSW ALRIKLYNKYV QGNAERLSVSSDDSCSEDGENPPRTHMNEGGYNTPVHH PVIKSENSVIKA GVRPESRASEDYVAPPKYKSSLLNRYLNDSLRHVMATNTT MMGKTRQRNH SGSFSSNEFEEELVSPSSSETEGNFVYRTDTLENGSNKG ESAANRKGPSK DDTYWKEINAAMALTNLAQGKDKLQGTTSCTIIQKSSHIAEV KTVKVPLVQ QF
Amino acids	1 to 352
Tags	GST tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab166113 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	ELISA Western blot
Form	Liquid
Additional notes	

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

May act as a morphogenetic regulator of cell adhesion.

Sequence similarities


Belongs to the TALE/IRO homeobox family.

Contains 1 homeobox DNA-binding domain.

Cellular localization

Nucleus.

Images



ab166113 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors