

Product datasheet

Recombinant Human MLLT11 protein (denatured) ab181903

1 Image

Description

Product name	Recombinant Human MLLT11 protein (denatured)
Purity	> 80 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>NP_006809</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMRDPVSS QYSSFLFWRM PIPELDLSEL EGLGLSDTAT YKVKDSSVGK MIGQATAADQ EKNPEGDGLL EYSTFNFWRA PIASHSFEL DLL
Predicted molecular weight	10 kDa
Amino acids	1 to 90
Tags	His tag N-Terminus
Description	Recombinant Human MLLT11 protein

Specifications

Our **Abpromise guarantee** covers the use of **ab181903** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	--


pH: 8.00

Constituents: 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 2.4% Urea

General Info

Tissue specificity	Highly expressed in thymus and all leukemic cell lines.
Involvement in disease	A chromosomal aberration involving MLLT11 is found in acute leukemias. Translocation t(1;11) (q21;q23) with KMT2A/MLL1.
Sequence similarities	Belongs to the MLLT11 family.

Images



15% SDS-PAGE analysis of 3ug ab181903 stained with Coomassie Blue.

SDS-PAGE - Recombinant Human MLLT11 protein
(denatured) (ab181903)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors