

Recombinant Human NAIF1 protein ab181935

1 Image

Description

Product name	Recombinant Human NAIF1 protein
Purity	> 90 % SDS-PAGE. ab181935 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<a href="#">Q69YI7</a>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSMAPPAKKRKMNFSE REVEIVEELELK KHLLVNHFNAGVPLAAKSAAWHGILRRVNAVATCRREL EVKKKWSDLKT EVRRKVAQVRAAVEGGEAPGPTTEEDGAGGPGTGGGSG GGGPAVAPVLLTP MQQRICNLLGEATIISLPSTTEIHPVALGPSATAAAATVTLTQ IPTETTY HTLEEGVVEYCTAEAPPPLPPETPVDMMQAHDTSVKPQ ALKSRIALNSA KLIQEQRVTNLHVKEIAQHLEQQNDLLQMIRRSQEVQACA QERQAQAMEG TQAALSVLIQVLRPMIKDFRRYLQSNTPAPASDPGQVA QNGQPDSIIQ
Predicted molecular weight	38 kDa including tags
Amino acids	1 to 327
Tags	His tag N-Terminus
Additional sequence information	NP_931045

Specifications

Our **Abpromise guarantee** covers the use of **ab181935** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Mass Spectrometry SDS-PAGE
<b>Mass spectrometry</b>	MALDI-TOF
<b>Form</b>	Liquid


## Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 30% Glycerol (glycerin, glycerine), 0.88% Sodium chloride, 0.32% Tris HCl, 0.02% DTT
------------------------------	--

## General Info

<b>Function</b>	Induces apoptosis.
<b>Tissue specificity</b>	Widely expressed.
<b>Sequence similarities</b>	Belongs to the NAIF1 family.
<b>Cellular localization</b>	Nucleus.

## Images



15% SDS-PAGE analysis of 3 µg ab181935.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors