

Product datasheet

Recombinant Human Nectin 2 protein (Tagged) (Biotin) ab271401

[1 Image](#)

Description

Product name	Recombinant Human Nectin 2 protein (Tagged) (Biotin)
Purity	>= 90 % SDS-PAGE.
Expression system	HEK 293 cells
Accession	<u>Q92692</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	QDVRVQVLPEVRGQLGGTVELPCHLLPPVPGLYISLVTW QRPDAPANHQN VAAFHPKMGPSFSPKPGSERLSFVSAKQSTGQDTEAE LQDATLALHGLT VEDEGNYTCEFATFPKGSVRGMTWLRVIAKPKNQAEQK VTFSQDPTTVA LCISKEGRPPARISWLSSLDWEAKETQVSGTLAGTVTVTS RFTLVPSGRA DGVTVTCKVEHESFEERALIPVTLVRYPPPEVSISGYDDN WYLGRTDNL SCDVRSNPEPTGYDWSSTSGTFPTSAVAQGSQQLVIHAVD SLFNNTTFVCTV TNAVGMGRAEQVIFVRETPNTAGAGATGG
Predicted molecular weight	38 kDa including tags
Amino acids	32 to 360
Tags	His tag C-Terminus
Conjugation	Biotin

Specifications

Our **Abpromise guarantee** covers the use of **ab271401** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

General Info


Function Probable cell adhesion protein.

Tissue specificity Ubiquitous.

Sequence similarities Belongs to the nectin family.
Contains 2 Ig-like C2-type (immunoglobulin-like) domains.
Contains 1 Ig-like V-type (immunoglobulin-like) domain.

Cellular localization Cell membrane.

Images



SDS-PAGE analysis of ab271401.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors