

## Product datasheet

# Recombinant Human NUR77 protein ab152448

[1 Image](#)

### Description

<b>Product name</b>	Recombinant Human NUR77 protein
<b>Expression system</b>	Wheat germ
<b>Accession</b>	<b><u>P22736</u></b>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	SFSPPTGPSPSLAQSPKLFP SQATHQLGEGESYSMP TAPGLAPTSPHL EGSGILDTPVTSTKARSGAPGGSEGRCAVCGDNASCQHY G
<b>Predicted molecular weight</b>	36 kDa including tags
<b>Amino acids</b>	191 to 280

### Specifications

Our **Abpromise guarantee** covers the use of **ab152448** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE Western blot ELISA
<b>Form</b>	Liquid
<b>Additional notes</b>	


### Preparation and Storage

<b>Stability and Storage</b>	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

## General Info

<b>Function</b>	Orphan nuclear receptor. May act concomitantly with NURR1 in regulating the expression of delayed-early genes during liver regeneration. Binds the NGF1-B response element (NBRE) 5'-AAAAGGTCA-3' (By similarity). May inhibit NF-kappa-B transactivation of IL2.
<b>Tissue specificity</b>	Fetal muscle and adult liver, brain and thyroid.
<b>Sequence similarities</b>	Belongs to the nuclear hormone receptor family. NR4 subfamily. Contains 1 nuclear receptor DNA-binding domain.
<b>Cellular localization</b>	Nucleus.

## Images



12.5% SDS-PAGE analysis of ab152448 stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors