

Recombinant Human PAQR8 protein ab164739

1 Image

Description

Product name	Recombinant Human PAQR8 protein
Expression system	Wheat germ
Accession	<u><b>Q8TEZ7</b></u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MTTAILERLSTLSVSGQRLRRLPKILEDGLPKMPCTVPETD VPQLFREPY IRTGYRPTGHEWRYFFSLFQKHNEVVNVWTHLLAALAVL LRFWAFAEAE ALPWASTHSLPLLLFILSSITYLTCSLLAHLLQSKSELSHYTF YFVDYVG VSVYQYGSALAHFFYSSDQAWYDRFWLFFLPAAAFCGW LSCAGCCYAKYC YRRPYPMRKICQVVPAGLAFILDISPVAHRVALCHLAGCQ EQAAWYHTL QILFFLVSAFFFCPVPEKYFPGSCDIVGHGHQIFHAFLSIC TLSQLEAI LLDYQGRQEIFLQRHGPLSVHMACLSFFFLAACSAATAAL LRHKVKARLT KKDS
Amino acids	1 to 354
Tags	GST tag N-Terminus

Specifications

Our <b>Abpromise guarantee</b> covers the use of <b>ab164739</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	ELISA  Western blot
Form	Liquid
Additional notes	This product was previously labelled as Progesterin Receptor Beta.

## Preparation and Storage

### Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info


### Relevance

MPRb (Membrane progesterin receptor beta) is a steroid membrane receptor and binds progesterone. MPRb may be involved in oocyte maturation. MPRb is an integral membrane protein and is expressed in brain and testis. It belongs to the ADIPOR family.

### Cellular localization

Cell Membrane. Multi-pass membrane protein. Note: Colocalizes with a lysosomal protein cathepsin D.

## Images



ab164739 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors