

Product datasheet

Recombinant human PARP11 protein ab79742

[3 Images](#)

Description

Product name	Recombinant human PARP11 protein
Biological activity	Specific Activity: 23 U/mg.
Purity	> 80 % SDS-PAGE. Affinity purified.
Expression system	Baculovirus infected Sf9 cells
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Specifications

Our **Abpromise guarantee** covers the use of **ab79742** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Functional Studies SDS-PAGE Western blot
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.0462% (R*,R*)-1,4-Dimercaptobutan-2,3-diol, 0.395% Tris HCl, 0.05% Tween, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------------------------	--

General Info

Relevance	Poly(ADP-ribosyl)ation is a DNA-damage-dependent post-translational modification of histones.
------------------	---

Poly(ADP-ribose) polymerases (PARPs) are a family of 18 proteins, encoded by different genes and displaying a conserved catalytic domain.

Images






Image showing specific activity of ab79742.

Functional Studies - Recombinant human PARP11 protein (ab79742)



10% SDS-PAGE showing ab79742 at approximately 65.6kDa (4 μ g).

SDS-PAGE - Recombinant human PARP11 protein (ab79742)



Western blot - Recombinant human PARP11 protein (ab79742)

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors