

Recombinant Human PARVA protein - BSA and Azide free ab173045

Description

Product name	Recombinant Human PARVA protein - BSA and Azide free
Purity	> 95 % SDS-PAGE. Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
Endotoxin level	< 1.000 Eu/μg
Expression system	Escherichia coli
Accession	<u>Q9NVD7</u>
Protein length	Full length protein
Animal free	No
Carrier free	Yes
Nature	Recombinant
Species	Human
Predicted molecular weight	43 kDa including tags
Amino acids	1 to 372
Tags	His tag C-Terminus
Description	Recombinant Human PARVA protein (BSA and azide free)

Specifications

Our **Abpromise guarantee** covers the use of **ab173045** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	HPLC SDS-PAGE
Form	Liquid
Additional notes	This product was previously labelled as Parvin alpha

Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	--

pH: 7.40

Constituents: 40% Glycerol (glycerin, glycerine), 0.88% Sodium chloride, 0.61% Tris

Supplied as a 0.2 µm Ø filtered solution.

General Info

Function	Probably plays a role in the regulation of cell adhesion and cytoskeleton organization. Plays a role in ciliogenesis.
Tissue specificity	Widely expressed, with highest levels in heart, skeletal muscle, kidney and liver.
Sequence similarities	Belongs to the parvin family. Contains 2 CH (calponin-homology) domains.
Cellular localization	Cell junction > focal adhesion. Cell membrane. Cytoplasm > cytoskeleton. Constituent of focal adhesions.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors