

# Recombinant Human Pellino 1 protein ab127100

## Description

<b>Product name</b>	Recombinant Human Pellino 1 protein
<b>Purity</b>	= 80 % SDS-PAGE. Purified via His tag
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<b><u>Q96FA3</u></b>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Predicted molecular weight</b>	29 kDa
<b>Amino acids</b>	145 to 404
<b>Tags</b>	His-DHFR tag N-Terminus

## Specifications

Our **Abpromise guarantee** covers the use of **ab127100** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE
<b>Form</b>	Lyophilized

## Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
<b>Reconstitution</b>	Reconstitute with water to desired concentration

## General Info

<b>Function</b>	Scaffold protein involved in the IL-1 signaling pathway via its interaction with the complex containing IRAK kinases and TRAF6. Required for NF-kappa-B activation and IL-8 gene expression in response to IL-1.
-----------------	--

## Sequence similarities

Belongs to the pellino family.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors