

Product datasheet

Recombinant Human Plunc protein ab139617

1 Image

Description

Product name	Recombinant Human Plunc protein
Purity	> 90 % SDS-PAGE. ab139617 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q9NP55</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MGSSHHHHHH SSGLVPRGSH MGSQFGGLPV PLDQTLPLNV NPALPLSPTG LAGSLTNALS NGLLSGGLLG ILENPLLDI LKPGGGTSGG LLGGLLGKVT SVIPGLNNII DIKVTDPQLL ELGLVQSPDG HRLYVTIPLG IKLQVNTPLV GASLLRLAVK LDITAEILAV RDKQERIHLV LGDCTHSPGS LQISLLDGLG PLPIQGLLDS LTGILNKVLP ELVQGNVCPL VNEVLRGLDI TLVHDMNML IHGLQFVIKV</p>
Predicted molecular weight	27 kDa including tags
Amino acids	20 to 256
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab139617** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.32% Tris HCl, 30% Glycerol, 1.17% Sodium chloride

General Info

Relevance

May be involved in the airway inflammatory response after exposure to irritants. May be associated with tumor progression. May play a role in innate immune responses of the upper airways.

Cellular localization

Secreted. Note: Found in the nasal mucus. Apical side of airway epithelial cells. Detected in nasal mucus

Images



15% SDS-PAGE analysis of ab139617 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors