

Product datasheet

Recombinant Human PNAD protein ab165200

1 Image

Description

Product name	Recombinant Human PNAD protein	
Expression system	Wheat germ	
Accession	<u>Q96AB6</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MPELLVEGRRVRLPQSAGDLVRAHPPLEERARLLRGQSVQ QVGPQGLLYVQ QRELAVTSPKDGSIILGSDDATTCMVL RHTGNGATCLT HCDGTDTKA EVPLIMNSIKSFSDDHAQCGRLEVHLVGGFSDDRQLSQKLT HQLLSEFDRQ EDDIHLVTL CVTELNDREENENHFPVIYGI AVNIKTAEYRAS FQDRGPE EQLRAARTLAGGPMISYDAETEQLRIGPYSWTFPHVDF WLHQDDKQIL ENLSTSP LAEPPHFVEHIRSTLMFLKKHPSPAHTLFSGNK ALLYKKNEDG LWEKISSPGS	
Amino acids	1 to 310	
Tags	GST tag N-Terminus	

Specifications

Our **Abpromise guarantee** covers the use of **ab165200** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid
Additional notes	This product was previously labelled as NTAN1.

Preparation and Storage


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Relevance	PNAD is required for the ubiquitin-dependent turnover of intracellular proteins that initiate with Met-Asn. This enzyme does not act on substrates with internal or C-terminal asparagines and does not act on glutamine residues in any position.
Cellular localization	Cytoplasmic

Images



ab165200 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors