

## Product datasheet

# Recombinant Human PREI3 protein ab101639

1 Image

### Description

---

<b>Product name</b>	Recombinant Human PREI3 protein
<b>Purity</b>	> 95 % SDS-PAGE. ab101639 was purified using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<u><b>Q9Y3A3</b></u>
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Predicted molecular weight</b>	28 kDa including tags
<b>Amino acids</b>	1 to 225
<b>Tags</b>	His tag N-Terminus

### Specifications

---

Our **Abpromise guarantee** covers the use of **ab101639** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE Mass Spectrometry
<b>Mass spectrometry</b>	MALDI-TOF
<b>Form</b>	Liquid

### Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 0.077% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 1.16% Sodium chloride
------------------------------	--

## General Info

---

<b>Sequence similarities</b>	Belongs to the MOB1/phocein family.
<b>Post-translational modifications</b>	Phosphorylated on serine residues.
<b>Cellular localization</b>	Cytoplasm > perinuclear region. Membrane. Golgi apparatus > Golgi stack membrane. In a perinuclear punctate pattern. Associated with membranes and the Golgi stacks.

---

## Images

---



ab101639 at 3 µg analysed by 15% SDS PAGE.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors