

Recombinant Human PRMT10 protein (Tagged) ab271722

[1 Image](#)

Description

Product name	Recombinant Human PRMT10 protein (Tagged)
Purity	= 90 % SDS-PAGE. Affinity purified.
Expression system	Baculovirus infected Sf9 cells
Accession	<u>Q6P2P2</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>SNSRPRSRR DAGGGAGAAG RDELVSRSLQ SAEHCLGVQD FGTAYAHYLL VLSLAPELKH DVKETFYQYL FRWAEELDAL SRIQDLLGCV EQALELFPDD EVICNSMGEH LFRMGFRDEA AGYFHKAVKL NPDFSDAKEN FYRVANWLVE RWHFIMLNDT KRNTIYNAAI QKAVCLGSKS VLDIGAGTGI LSMFAKKAGA HSVYACELSK TMYELACDVV AANKMEAGIK LLHTKSLDIE IPKHIPERVS LVTETVDAG LFGEVIVESL IHAWHEHLLQ PKTKGESANC EKYGKVIPAS AVIFGMAVEC AEIRRHHRVG IKDIAGIHLN TNVKFQSPAY SSVDTTEETIE PYTTEKMSRV PGGYLALTEC FEIMTVDFNN LQELKSLATK KPDKIGIPVI KEGILDAMV WFLVQLDDEH SLSTSPSEET CWEQAVYPVQ DLADYWIKPG DHVMMEVSCQ DCYLRIQISIS VLGLECEMDV AKSFTQNKDL LSLGNEAELC SALANLQTSK PDAVEQTCIL ESTEIALLN IPYHEGFKMA MSKVLSSLTP EKLYQTMETH CQNEMSSGTG QSNTVQNILE PFYVLDVSEGF FSVLPVIAGT LGQVKPYSSV EKDQHRIALD LISEANHFPE ETLEFWLRHV EDESAMLQRP KSDKLWSIII LDVIEPSGLI QQEIMEKAAI SRCLLQSGGK IFPQVYLMFG LLVESQTLLE ENAVQGTERT LGLNIAPFIN QFQVPIRVFL DLSSLPCIPL SKPVELLRLD LMTPLYNLSN REVKVYVCKS GRLTAIPFWY HMYLDEEIRL DTSSEASHWK</p>

Predicted molecular weight	95 kDa including tags
Amino acids	2 to 845
Tags	His tag C-Terminus , DDDDK tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab271722** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.63% Tris HCl, 0.64% Sodium chloride, 0.02% Potassium chloride, 0.04% Tween, 20% Glycerol (glycerin, glycerine), 0.06% TCEP 100ug/ml flag peptide
------------------------------	---

General Info

Sequence similarities	Belongs to the protein arginine N-methyltransferase family. Contains 3 TPR repeats.
------------------------------	--

Images



SDS-PAGE analysis of ab271722.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors