

Recombinant Human PSG2 protein ab185849

Description

Product name	Recombinant Human PSG2 protein
Purity	> 95 % SDS-PAGE. Assessed also by SEC-HPLC
Endotoxin level	< 1.000 Eu/μg
Expression system	Mammalian
Accession	<u>P11465</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	QVTIEAQPPKVSEGKDVLLL VHNLPQNL TGYWYKGQIRDL YHYITSYVV DGQIIYGPA YSGRE TAYSNASLLIQNVTREDAGSYTLHIIKRGDGT RGV TGYFTFT LYLETPKP SSSSNLNPREAMET VILTCDPETPDTSYQWWMNGQSLPMTHRFQLSE TNRTLFL FGVTKYTAGP YECEIRNSGSASRSD PVTLNLLHGPDLPRIHPSYTNYRSGDNL YLSCFANSNPPA QYSWTINGKF QQSGQNLFI PQITK HSGLYVCSVRNSATGEESSTSLTVKVSASTRIGLLPLL NPT VDHHHHHH
Predicted molecular weight	35 kDa including tags
Amino acids	35 to 335
Tags	His tag C-Terminus
Additional sequence information	Mature form.

Specifications

Our **Abpromise guarantee** covers the use of **ab185849** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	HPLC SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 99% Phosphate Buffer, 0.88% Sodium chloride

Supplied as a 0.2 µM filtered solution

General Info

Sequence similarities

Belongs to the immunoglobulin superfamily. CEA family.

Contains 2 Ig-like C2-type (immunoglobulin-like) domains.

Contains 1 Ig-like V-type (immunoglobulin-like) domain.

Cellular localization

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors