

Recombinant Human RASAL1 protein ab160027

1 Image

Description

Product name	Recombinant Human RASAL1 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MAKSSSLNVRVVEGRALPAKDVSGSSDPYCLVKVDDEV VARTATVWRSLG PFWGEEYTVHLPDLDFHQLAFYVLDEDTVGHDDIIGKISLSR EAITADPRG IDSWINLSRVPDAEVQGEICLSVQMLEDGQGRCLRCHVL QARDLAPRDI SGTSDPFARVFWGSQSLETSTIKKTRFPHWDEVLELREM PGAPSPLRVEL WDWDMVGKNDFLGMVEFSPKTLQQKPPKGWFRLLPFP RAEEDSGGNLGAL RVKVRLIEDRVLP SQCYQPLMELLMESVQGPAEEDTASP LALLEELTLGD CRQDLATKLVKLFLGRGLAGHFLDYLTRREVARTMDPNTL FRSNSLASKS MEQFMKLVGMPYLHEVLKPVISRVFEEKKYMELDPCKMD LGRTRRISFKG ALSEEQMRETSLGLLTGYLGPIVDANGSVGRCPAMRLA FKQLHRRVEE RFPQAEHQDVKYLAISGFLFLRFFAPAILTPKLFDLRDQHA DPQTSRSL LLAKAVQSIGNLGQQLGQKELWMAPLHPFLLQCVSRVR DFLDRLVDVDG DEAGVPARALFPPSAMREGYLLKRKEEPAGLATRFAFKK RYVWLSGETL SFSKSPEWQVVTQDGTGALHTTYLQCKNVNELNQWLSAL RKASAPNPKNL AACHPGAFRSARWTCCLQAERSAAGCSRTHSAVTLGDW SDPLDPDAEAQT VYRQLLLGRDQLRLKLLSDNMDTTLEADTGACPEVLAR

Amino acids	1 to 776
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab160027** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot
	ELISA
Form	Liquid

Additional notes


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Function	Probable inhibitory regulator of the Ras-cyclic AMP pathway.
Tissue specificity	Highly expressed in thyroid and adrenal medulla, lower expression in brain, spinal cord and trachea.
Sequence similarities	Contains 1 Btk-type zinc finger. Contains 2 C2 domains. Contains 1 PH domain. Contains 1 Ras-GAP domain.

Images



ab160027 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors