

Product datasheet

Recombinant Human Secretogranin 3 protein ab167866

1 Image

Description

Product name	Recombinant Human Secretogranin 3 protein
Purity	> 85 % SDS-PAGE. ab167866 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q8WXD2</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MGSSHHHHHH SSGLVPRGSH MGSFPPKPGGS QDKSLHNREL SAERPLNEQIAEAEEDKIKK TYPENKPGQ SNYSFVDNLN LLKAITEKEK IEKERQSIRS SPLDNKLNVE DVDSTKNRKL IDDYDSTKSG LDHKFQDDPD GLHQLDGTPL TAEDIVHKIA ARIYEENDRA VFDKIVSKLL NLGLITESQA HTLEDEVAEV LQKLISKEAN NYEEDPNKPT SWTENQAGKI PEKVTPMAAI QDGLAKGEND ETVSNTLTLT NGLERRTKTY SEDNFEELQY FPNFYALLKS IDSEKEAKEK ETLITIMKTL IDFKMMVKY GTISPEEGVS YLENLDEMIA LQTKNKLEKN ATDNISLFP APSEKSHEET DSTKEEAAKM EKEYGSLKDS TKDDNSNPGG KTDEPKGKTE AYLEAIRKNI EWLKKHDKKG NKEDYDLSKM RDFINKQADA YVEKGILDKE EAEAIKRIYS SL</p>
Predicted molecular weight	53 kDa including tags
Amino acids	20 to 468
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab167866** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.32% Tris-HCl buffer, 0.02% DTT, 20% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

General Info

Tissue specificity Expressed in brain, heart, kidney, liver and skeletal muscle.

Cellular localization Secreted. Neuroendocrine and endocrine secretory granules.

Images



15% SDS-PAGE analysis of 3µg ab167866.

SDS-PAGE - Recombinant Human Secretogranin 3 protein (ab167866)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors