

Recombinant Human SERBP1 protein ab153267

1 Image

Description

Product name	Recombinant Human SERBP1 protein		
Expression system	Wheat germ		
Protein length	Full length protein		
Animal free	No		
Nature	Recombinant		
Species	Human		
Sequence	MPGHLQEGFGCVVTNRFDQLFDDSDPFVFLKAAENKK KEAGGGGVGGPG AKSAAQAAAQTSNAAGQLRKESQKDRKNPLPPSVGV VDKKEETQPPVA LKKEGIRRVGRRPDQQLQGEGKIIDRRPERRPPRERRFEK PLEEKGEGGE FSVDRPIIDRPIRGRGGLGRGRGGRGRGMGRGDGFDSRG KREFDRHSGSD RSGLKHEDKRGSGSHNWGTVKDELTESPKYIQKQISYNY SDLDQSNVTE ETPEGEEHHPVADTENKENEVEEVKEEGPKEMTLDEWK AIQNKDRAKVEF NIRKPNEGADGQWKKGFVLHKSKEEAHAEDSVMDHHF RKPANDITSQLE INFGDLGRPGRGGRGGRGGCGRGRPNRGSRTDKSSAS APDVDDPEAFPA LA		
Amino acids	1 to 402		
Tags	GST tag N-Terminus		

Specifications

Our **Abpromise guarantee** covers the use of **ab153267** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot
Form	Liquid

Additional notes

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

May play a role in the regulation of mRNA stability. Binds to the 3'-most 134 nt of the SERPINE1/PAI1 mRNA, a region which confers cyclic nucleotide regulation of message decay.


Tissue specificity

Expressed at high level in the heart, skeletal muscle and kidney, and at low levels in placenta, liver and brain.

Cellular localization

Cytoplasm. Nucleus. Cytoplasm > perinuclear region. Also found in perinuclear regions.

Images



SDS-PAGE - Recombinant Human SERBP1 protein
(ab153267)

ab153267 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors