

Product datasheet

Recombinant Human SHK protein ab161754

1 Image

Description

Product name	Recombinant Human SHK protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MAARPITLGIDLGTTSVKAALLRAAPDDPSGFAVLASCAR AARAEAAVES AVAGPQGREQDVSRILQALHECLAALPRPQLRSVVGIGVS GQMHGVVFWK TGQGCEWTEGGITPVFEPRAVSHLVTWQDGRCSSEFLAS LPQPKSHLSVA TGFGCATIFWLLKYRPEFLKSYDAAGTIHDYVAMLCGLPR PLMSDQNA SWGYNFNTQSQSWNVETLRSSGFPVHLLPDIAEPGSVAGR TSHMWFEIPKG TQVGVALGDLQASVYSCMAQRDVAVLNISTSVQLAASMP SGFQPAQTPDP TAPVAYFPYFNRTYLGVAASLNNGNVLATFVHMLVQWMA DLGLEVEESTV YSRMIQAAVQQRDTHLTITPTVLGERHLPDQLASVTRISS DLSLGHVTR ALCRGIVQNLHSMLPIQQLQEWGVERVMGSGSALS RNDV LKQEVQRAFP PMSFGQDVDAAVGAALVMLRRHLNQKES
Amino acids	1 to 478
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab161754** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Western blot

	ELISA
Form	Liquid
Additional notes	This product was previously labelled as CARKL.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Tissue specificity	Strongly expressed in liver, kidney and pancreas. Expressed at lower levels in placenta and heart. Very weakly expressed in lung and brain.
Sequence similarities	Belongs to the FGGY kinase family.
Cellular localization	Cytoplasm.

Images



ab161754 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors