

Product datasheet

Recombinant Human SLFN11 protein ab164865

[1 Image](#)

Description

Product name	Recombinant Human SLFN11 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Sequence

MEANQCPLVVEPSYPDLVINVGEVTLGEEENRKKLQKIQRD
QEKERVMAAA
CALLNSGGGVIRMAKKVEHPVEMGLDLEQSLRELIQSSDL
QAFFETKQQG
RCFYIFVKSWSGGPFPEDRSVKPRCLSLSSSLYRRSETSV
RSMDSREAFCA
FLKTKRKPKILEEGPFHKKHKGVYQELPNSDPADPNSDPA
DLIFQKDYLE
YGEILPFPESQLVEFKQFSTKHFQEYVKRTIPEYVPAFANT
GGGYLFIGV
DDKSREVLGCAKENVDPDSLRRKIEQAYKLPCVHFCQP
QRPITFTLKIV
DVLKRGELYGYACMIRVNPFCFAVFSEAPNSWIVEDKYV
CSLTTEKWVGM
MTDTPDLLQLSEDFECQLSLSSGPPLSRPVYSKKGLEH
KKELQQLLFSV
PPGYLRYTPESLWRDLISEHRGLEELINKMQPFFRGILIFS
RSWAVDLN
LQEKPGVICDALLIAQNSTPILYTLREQDAEGQDYCTRATF
TLKQKLVN
MGGYTGKVCVRAKVLCLSPESAEAEAAVSPMDYPASY
SLAGTQHMEAL
LQSLVIMLLGFRSLLSDQLGCEVLNLLTAQQYEIFSRSLRK
NRELFVHGL
PGSGKTIMAMKIMEKIRNVFHCEAHRILYVCENQPLRNFISD
RNICRAET
RKTFLRENFEBHIQHVIDEAQNFRTEDGDWYGKAKSITRRA
KGGPGLWI
FLDYFQTSHLDCSGLPPLSDQYPREELTRIVRNADPIAKYL

QKEMQVIRS
NPSFNIPTGCLEVFPEAEWSQGVQGTLRICKKYLTV EQIMTC
VADTCRRFF
DRGYSPKDVAVLVSTAKEVEHYKYELLKAMRKKRVVQLS
DACDMLGDHV
LDSVRRFSGLERSIVFGIHPRTADPAILPNVLICLASRAKQH
LYIFPWGG H

Amino acids 1 to 901
Tags GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab164865** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Western blot
ELISA
Form Liquid

Additional notes


Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Sequence similarities Belongs to the Schlafen family.

Images



ab164865 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human SLFN11 protein
(ab164865)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors