

Product datasheet

Recombinant Human SPAG7 protein ab168077

1 Image

Description

Product name	Recombinant Human SPAG7 protein
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u><b>O75391</b></u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMADLLGS ILSSMEKPPS LGDQETRRKA REQAARLKKL QEQEKKQKVE FRKRMEKEVS DFIQDSGQIK KKFQPMNKIE RSILHDVVEV AGLTSFSFGE DDDCRYVMIF KKEFAPSDEE LDSYRRGEEW DPQKAEERKR LKELAQRQEE EAAQQGPVVV SPASDYKDKY SHLIGKGAAC DAAHMLQANK TYGCVPVANK RDTRSIEEAM NEIRAKKRLR QSGEELPPTS
Predicted molecular weight	28 kDa including tags
Amino acids	1 to 227
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab168077** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

## Preparation and Storage

### Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.50

Constituents: 0.02% DTT, 0.32% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

## General Info

### Tissue specificity

Detected in fetal brain.


### Sequence similarities

Contains 1 R3H domain.

### Cellular localization

Nucleus.

## Images



3 µg ab168077 on 15% SDS-PAGE

SDS-PAGE - Recombinant Human SPAG7 protein  
(ab168077)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors