

Product datasheet

Recombinant Human SRSF3 protein ab114675

[1 References](#) [1 Image](#)

Description

Product name	Recombinant Human SRSF3 protein
Expression system	Wheat germ
Accession	<u>P84103</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MHRDSCPLDCKVYVGNLGNNGNKTELERAFGYGPLRSV WVARNPPGFAF VEFEDPRDAADAVRELDGRTLCGCRVRVELSNGEK
Predicted molecular weight	35 kDa including tags
Amino acids	1 to 85

Specifications

Our **Abpromise guarantee** covers the use of **ab114675** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot SDS-PAGE ELISA
Form	Liquid
Additional notes	This product was previously labelled as SFRS3.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------------------------	---

General Info

Function	May be involved in RNA processing in relation with cellular proliferation and/or maturation.
Sequence similarities	Belongs to the splicing factor SR family. Contains 1 RRM (RNA recognition motif) domain.
Post-translational modifications	Extensively phosphorylated on serine residues in the RS domain.
Cellular localization	Nucleus.

Images



SDS-PAGE analysis of ab114675 on a 12.5% gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors