

Recombinant Human TBC1D22B protein ab202146

1 Image

Description

Product name	Recombinant Human TBC1D22B protein
Purity	> 80 % SDS-PAGE. ab202146 was purified by conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q9NU19</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAAENSK QFWKRSAKLP GSIQPVYGAQ HPPLDPRLTK NFIKERSKVN TVPLKNKKAS SFHEFARNTS DAWDIGDDEE EDFSSPSFQT LNSKVALATA AQVLENHSKL RVKPERSONQST TSDVPANYKV IKSSSDAQLS RNSSDTCLRN PLHKQQSLPL RPIIPLVARI SDQNASGAPP MTVREKTRLE KFRQLLSSQN TDLDELKCS WPGVPREVRP ITWRLLSGYL PANTERRKLT LQRKREEYFG FIEQYDSRN EEHHQDTYRQ IHIDIPRTNP LIPLFQQPLV QEIFERILFI WAIRHPASGY VQGINDLVTP FFVVFLSEYV EEDVENFDVT NLSQDMLRSI EADSFWCMSK LLDGIQDNYT FAQPGIQKKV KALEELVSRI DEQVHNHFRR YEVEYLQFAF RWMNNLLMRE LPLRCTIRLW DTYQSEPEGF SHFHLVCAA FLIKWRKEIL DEEDFQGLLM LLQNLPTIHW GNEEIGLLLA EAYRLKYMFA DAPNHYYR
Predicted molecular weight	62 kDa including tags
Amino acids	1 to 505
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab202146** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 7.40

Constituents: 30% Glycerol (glycerin, glycerine), 0.02% DTT, 69% PBS

General Info

Relevance TBC1D22B may act as a GTPase activating protein for Rab family protein(s).

Images



3ug by SDS-PAGE under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human TBC1D22B protein (ab202146)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors