

Product datasheet

Recombinant Human TDGF1P3 protein (His tag) ab267981

Description

Product name	Recombinant Human TDGF1P3 protein (His tag)
Purity	> 85 % SDS-PAGE.
Expression system	Mammalian
Accession	<u>P51864</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	LGHQEFARPSRGDLAFRDDSIWPQEEPAIRPRSSQRVLP MGIQHSKELNR TCCLNGGTCMLSEFCACPPSFYGRNCEHDVRKENCGSV PHDTWLPPKKCSL CKCWHGQLRCFPQAFLPGCDGLVMDEHLVASRTPELPP SARTTTFMLAGICLSIQSY
Amino acids	31 to 188
Tags	His tag C-Terminus
Additional sequence information	Mature chain

Specifications

Our **Abpromise guarantee** covers the use of **ab267981** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	----------------------------------------------------------------------------------------------------------------------------------

pH: 7.2
Constituents: PBS, 6% Trehalose

Lyophilized from.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors