

Product datasheet

Recombinant Human TNFAIP2 protein ab159728

1 Image

Description	
Product name	Recombinant Human TNFAIP2 protein
Expression system	Wheat germ
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	YILANADTIQHFCTQHGSPATWLQPALPTLAEIIRLQDPSAIK IEVATYA TCYPDFSKGHL SAILAIKGNLSNSEVKRIRSILDVSMGAQE PSRPLFSL
Amino acids	552 to 650
Tags	GST tag N-Terminus


Specifications	
Our Abpromise guarantee covers the use of ab159728 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	ELISA Western blot
Form	Liquid
Additional notes	

Preparation and Storage	
Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl

General info

Function	May play a role as a mediator of inflammation and angiogenesis.
Sequence similarities	Belongs to the SEC6 family.
Developmental stage	Differentially expressed in development and capillary tube-like formation in vitro.

Images



ab159728 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors