

Product datasheet

Recombinant Human TPD52L1 protein ab48771

1 Image

Description	
Product name	Recombinant Human TPD52L1 protein
Purity	> 90 % SDS-PAGE. TPD52L1 (1-128aa) with N-terminal His-tag was overexpressed in E. coli and purified by using conventional chromatography techniques.
Expression system	Escherichia coli
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MEAQAQGLLE TEPLQGTDED AVASADFSSM LSEEEKEELK AELVQLEDEI TTLRQVLSAK ERHLVEIKQK LGMNLMNELK QNFSKSWHDM QTTTAYKKTH ETLSHAGQKA TAAFSNVGTA ISKKFGDMRR K
Amino acids	1 to 128

Specifications	
Our <b>Abpromise guarantee</b> covers the use of <b>ab48771</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE
Form	Liquid


Preparation and Storage	
Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.40 Constituent: PBS

General Info	
--------------	--

## Sequence similarities

Belongs to the TPD52 family.

## Images



14% SDS-PAGE gel loaded with ab48771, recombinant human TPD52L1 fusion protein (his-tagged). Predicted molecular weight 16.8 kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors