

Recombinant Human TRIP-Br2 protein ab180305

1 Image

Description	
Product name	Recombinant Human TRIP-Br2 protein
Purity	> 90 % SDS-PAGE. ab180305 was purified by using conventional chromatography techniques
Expression system	Escherichia coli
Accession	<u>Q14140</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSM LGKGGKRKFDEHE DGLEGKIVSPCDG PSKVSYTLQRQTIFNISLMKLYNHRPLTEPSLQKTVLINNML RRIQEELK QEGSLRPMFTPSSQPTTEPSDSYREAPPAFSHLASPSSH PCDLGSTTPLE ACLT PASLLEDDDDTFCTSQAMQPTAPTKLSPPALLPEK DSFSSALDEIE ELCPTSTSTEATAATDSVKGTSSEAGTQKLDGPQESRA DDSKLMDSLPG NFEITTSTGFLTDLTLDDILFADIDTSMYDFDPCTSSSGTAS KMAPVSAD DLLKTLAPYSSQPVTPSQPFKMDLTELDHIMEVLVGS
Predicted molecular weight	36 kDa including tags
Amino acids	1 to 314
Tags	His tag N-Terminus
Additional sequence information	NP_055570

Specifications	
Our Abpromise guarantee covers the use of ab180305 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Mass Spectrometry


	SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid
Additional notes	Previously labelled as SERTAD2.

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.32% Tris HCl, 0.88% Sodium chloride, 10% Glycerol (glycerin, glycerine), 0.02% DTT
------------------------------	--

General Info

Images



15% SDS-PAGE analysis of ab180305 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors