

Recombinant Human TXNL1 protein ab103050

1 Image

Description

Product name	Recombinant Human TXNL1 protein
Purity	> 95 % SDS-PAGE. ab103050 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>O43396</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSH MGVVKPVGSDPDFQPEL SGAGSRLAVVKFT MRGCGPCLRIAPAFSSMSNKYPQAVFLEVDVHQCQGTAATNNISATPTFL FFRNKVRIDQYQGADAVGLEEKIKQHLENDPGSNEDTDIP KGYMDLMPFI NKAGCECLNESDEHGFDNCLRKDTTFLESDCDEQLLITV AFNQPVKLYSM KFQGPDNGQGPKYVKIFINLPRSMDFEEAERSEPTQALETEDDIKEDGI VPLRYVKFQNVNSVTIFVQSNQGEEETTRISYFTFIGTPVQATNMNDFKR VVGKKGESH
Predicted molecular weight	34 kDa including tags
Amino acids	1 to 289
Tags	His tag N-Terminus

Specifications

Our Abpromise guarantee covers the use of ab103050 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF

Form Liquid


Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
pH: 8.00
Constituents: 0.077% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 1.16% Sodium chloride

General Info

Tissue specificity Ubiquitous.
Sequence similarities Contains 1 PITH domain.
Contains 1 thioredoxin domain.
Cellular localization Cytoplasm.

Images



SDS-PAGE analysis of ab103050

SDS-PAGE - Recombinant Human TXNL1 protein
(ab103050)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors