

Recombinant Human Ube2L6 protein ab173041

Description

Product name	Recombinant Human Ube2L6 protein
Purity	> 95 % SDS-PAGE. Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
Endotoxin level	= 1.000 Eu/μg
Expression system	Escherichia coli
Accession	<u>O14933</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MMASMRVVKELEDLQKKPPPYLRNLSSDDANVLVWHAL LLPDQPPYHLKA FNLRFSPPEYFPKP PMIKFTTKMYHPNVDENGQICLPISSENWKPCCTKTCQVLEA LNVLVNRP NIREPLRMDLADLLT QNPELFRKNAAEEFTLRFGVDRPSLEHHHHHH
Predicted molecular weight	18 kDa
Amino acids	1 to 153
Tags	His tag C-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab173041** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	HPLC SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00
-----------------------	--

Constituents: 1.19% HEPES, 0.58% Sodium chloride

Supplied as a 0.2 µM filtered solution.

General Info

Function	Catalyzes the covalent attachment of ubiquitin or ISG15 to other proteins. Functions in the E6/E6-AP-induced ubiquitination of p53/TP53.
Tissue specificity	Present in natural killer cells (at protein level).
Pathway	Protein modification; protein ubiquitination.
Sequence similarities	Belongs to the ubiquitin-conjugating enzyme family.
Post-translational modifications	ISGylated.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors