

Product datasheet

Recombinant Human USP35 protein ab127080

Description

| | |
|-----------------------------------|--|
| Product name | Recombinant Human USP35 protein |
| Purity | > 90 % SDS-PAGE. Purified via His tag |
| Expression system | Escherichia coli |
| Accession | <u>A2RRA6</u> |
| Protein length | Protein fragment |
| Animal free | No |
| Nature | Recombinant |
| Species | Human |
| Predicted molecular weight | 24 kDa |
| Amino acids | 56 to 264 |

Specifications

Our **Abpromise guarantee** covers the use of **ab127080** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|---------------------|-------------|
| Applications | SDS-PAGE |
| Form | Lyophilized |

Preparation and Storage

| | |
|------------------------------|--|
| Stability and Storage | Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride |
| Reconstitution | Reconstitute with water to desired concentration. |

General Info

| | |
|------------------------------|--|
| Tissue specificity | Expressed in testis, pancreas and skeletal muscle. |
| Sequence similarities | Belongs to the peptidase C19 family. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors