

Product datasheet

Recombinant Human VISTA protein (Tagged)
ab271380

1 Image

Description	
Product name	Recombinant Human VISTA protein (Tagged)
Purity	>= 90 % SDS-PAGE.
Expression system	HEK 293 cells
Accession	<u>Q9H7M9</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	FKVATPYSLYVCPEGQNVTL TCRLLGPVDK GHDVTFYKTW YRSSRGEVQT CSERRPIRNL TFQDLHLHHG GHQAANTSHD LAQRHGLESA SDHHGNFSIT MRNLTLLDSG LYCCLVVEIR HHHSEHRVHG AMELQVQTGK DAPSNCVVYP SSSQDSENT AA
Predicted molecular weight	21 kDa including tags
Amino acids	33 to 194
Tags	His tag C-Terminus , Avi tag C-Terminus
Additional sequence information	Extracellular domain.

Specifications	
Our Abpromise guarantee covers the use of ab271380 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	SDS-PAGE
Form	Liquid

Preparation and Storage	
-------------------------	--

Stability and Storage

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

General Info


Sequence similarities

Contains 1 Ig-like (immunoglobulin-like) domain.

Cellular localization

Membrane.

Images



SDS-PAGE analysis of 4 µg ab271380.

This protein runs at a higher MW due to glycosylation.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors