

Product datasheet

Recombinant Human ZBTB46 protein ab165467

1 Image

Description

Product name	Recombinant Human ZBTB46 protein
Expression system	Wheat germ
Accession	<u>Q86UZ6</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MNNRKEDMEIASHYRHLLRELNEQRQHGVLCDCVVEG KVFKAHKNVLL GSSRYFKTLYCQVQKTSEQATVTHLDIVTAQGFKAIIDFMY SAHLALTSR NVIEVMSAASFLQMTDIVQACHDFIKAALDISIKSDASDELA EFEIGASS SSSTEALISAVMAGRSISPWLARRTSPANSSGDSAIASCH DGGSSYGKED QEPKADGPDDVSSQPLWPGDVGYGPLRIKEEQVSPSQY GGSELPSAKDGA VQNSFSEQSAGDAWQPTGRRKNRKNKTVRHITQQVED DSRASSPVPSFL PTSGWPFSSRDSNADLSVTEASSSDSRGERAELYAQVE ELLGGEASYLG PPLTPEKDDALHQATAVANLRAALMSKNSLLSLKADVLG DDGSLLFEYLP RGAHSLSLNEFTVIRKKFKCPYCSFSAMHQCILKRHMRSR TGERPYPCEI CGKKFTRREHMKRHTLVHSDKKYVCKVCSRVFMSAAS VGIRHGSRRHGV CTDCAGRGMAGPLDHGGGGGEGSPEALFPGDGPYLE PEDPRGEAEELGE DDEGLAPEDALLADDKDEEDSPRPRSPPGGPKDFAWL S</p>
Amino acids	1 to 589
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab165467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications ELISA
Western blot

Form Liquid

Additional notes

Preparation and Storage


Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Relevance ZBTB46 may be involved in transcriptional regulation. It contains one BTB (POZ) domain and two C2H2 type zinc fingers.

Cellular localization Nucleus

Images



SDS-PAGE - Recombinant Human ZBTB46 protein
(ab165467)

ab165467 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors