

Recombinant Human ZNF133 protein ab127455

Description

Product name	Recombinant Human ZNF133 protein
Purity	> 95 % SDS-PAGE. Purified via His tag
Expression system	Escherichia coli
Accession	<u>P52736-2</u>
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Predicted molecular weight	28 kDa
Amino acids	388 to 635
Tags	His-DHFR tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab127455** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
Reconstitution	Reconstitute with water to desired concentration.

General Info

Function	May be involved in transcriptional regulation as a repressor.
Tissue specificity	Seems ubiquitous. Seen in the heart, brain, placenta, lung, liver, skeletal muscle, kidney and

pancreas.

Sequence similarities

Belongs to the krueppel C2H2-type zinc-finger protein family.
Contains 15 C2H2-type zinc fingers.
Contains 1 KRAB domain.

Cellular localization

Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors