

Product datasheet

Recombinant Human ZNF581 protein ab169912

Description

Product name	Recombinant Human ZNF581 protein
Purity	> 90 % SDS-PAGE. ab169912 was expressed in E.coli as inclusion bodies. The final product was refolded and chromatographically purified.
Expression system	Escherichia coli
Accession	<u>Q9P0T4</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MASMTGGQQMGRGEFMLVLPSPCPQPLAFSSVETMEGP PRRTCRSPEPGP SSSIGSPQASSPPR PNHYLLIDTQGVPTVTLVDEESQREPGASGAPGQ KKCYSCPVCSRFEYMSYLQRHSITHSEVK PFECDICGKAFKRASHLA RHHSIHLAGGGRPHGCPLCPRRFRDAGELAQHSRVHSGE RPFQCPH CP RRFMEQNTLQKHTRWKHP
Predicted molecular weight	24 kDa including tags
Amino acids	1 to 197
Tags	T7 tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab169912** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00

Constituent: 0.32% Tris HCl

Contains NaCl, KCl, EDTA, arginine, DTT and glycerol.

General Info

Function	May be involved in transcriptional regulation.
Sequence similarities	Contains 4 C2H2-type zinc fingers.
Cellular localization	Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors