

Product datasheet

Recombinant Hydrophobin 2 protein (Tagged) ab236203

1 Image

Description

Product name	Recombinant Hydrophobin 2 protein (Tagged)	
Purity	> 85 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	<u>P79073</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Sequence	AVCPTGLFSNPLCCATNVLDLIGVDCKTPTIAVDTGAI FQA HCASKGSKP LCCVAPVADQALLCQKAIGTF	
Predicted molecular weight	13 kDa including tags	
Amino acids	16 to 86	
Tags	His tag N-Terminus	
Additional sequence information	Hypocrea jecorina (Trichoderma reesei). This product is the mature full length protein from aa 16 to 86. The signal peptide is not included. N-terminal 10xHis-B2M-JD-tag.	

Specifications

Our **Abpromise guarantee** covers the use of **ab236203** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.2 Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)
------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis with 5% enrichment gel and 15% separation gel of ab236203.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors