

Product datasheet

Recombinant Influenza A Hemagglutinin protein (His tag) (His tag) Strain Egypt/N03072/2010 ab217653

1 Image

Description

Product name	Recombinant Influenza A Hemagglutinin protein (His tag) (His tag) Strain Egypt/N03072/2010	
Purity	> 95 % SDS-PAGE.	
Expression system	HEK 293 cells	
Accession	<u>D5XIT5</u>	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Sequence	<p>QIC IGYHANNSTE QVDTIMEKNV TVTHAQDILE KTHNGKLCNL DGVKPLILRD CSVAGWLLGN PMCDEFLNVP EWSYVEKIN PANDLCYPGN FNDYEELKHL LSRINHFEDI QIIPKSSWSD HEASGVSSXC PYQGRSSFFR NVVWLTKKDN AYPTIKRSYN NTNQEDLLVL WGIHHPNDAA EQTRLYQNPT TYISVGTSTL NQRLVPKIAT RSKVNGQSGR MEFFWTILKS NDAINFESNG NFIAPENAYK IVKKGDSTIM KSELEYGN CN TKCQTPIGAI NSSMPFHNIH PLTIGECPKY VKSNRLVLAT GLRNSPQGER RRKKRGLFGA IAGFIEGGWQ GMVDGWYGYH HSNEQGSYGA ADKESTQKAI DGVTNKVN SIIDKMNTQFEA VGREFNLER RIENLNKKME DGFLDVWTYN AELLVLMENE RTLDFHDSNV KNLYDKVRLQ LRDNAKELGN GCFEFYHRCD NECMESVRNG TYDYPQYSEE ARLKREEISG VKLESIGTYQ</p>	
Amino acids	18 to 530	
Tags	His tag C-Terminus	
Additional sequence information	Influenza A (H5N1)(A/Egypt/N03072/2010) . Genebank accession: ADG21447.	
Description	Recombinant Influenza A Hemagglutinin protein (His tag)	

Specifications

Our **Abpromise guarantee** covers the use of **ab217653** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE


Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.

Constituent: 100% PBS

Images



SDS-PAGE analysis of ab217653.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors